Development of a methodology to enable non-linear in vitro-in vivo correlation for complex long-acting injections \

Jake Dickinson, Claire Patterson, Parmesh Gajjar and Paul A Dickinson

Seda Pharmaceutical Development Services® , Alderley Edge, SK10 4TG, UK



- Linear in vitro in vivo correlation (IVIVC) methods for both oral and non-oral dosage forms are reported [1-5], however linear IVIVC may be inappropriate for complex parenterals
- Long-acting injections (LAIs), with complex release profiles, likely require a **non-linear** relationship to correlate accelerated *in vitro* release to real time *in vivo* release

Methods

To demonstrate the steps required for non-linear IVIVC, we simulated datasets for accelerated in vitro dissolution and PK profiles for three different formulations that are typical of a parenteral PLGA microsphere product

Steps to perform non-linear IVIVC

- 1. Deconvolution of the absorption profile of each formulation
- 2. Model dissolution profiles and calculate scaled in vitro timepoints
- 3. Creating a Levy plot, by plotting in vivo timepoints against the scaled in vitro timepoints. The IVIVC is said to be non-linear if the Levy plot is best described by a non-linear function
- 4. Generation of a scaled dissolution profile
- Simulation of PK profiles using scaled dissolution. Assessment of the resulting IVIVC against the guideline criteria set by the FDA [1]





Figure 1 Flow chart of IVIVC methodology

Results

- Linear and non-linear IVIVC performed on 3 formulations. Example of IVIVC steps on one formulation is shown in Figures 2 & 3
- For the non-linear approach, high-order polynomials were used to describe the in vitro-in vivo relationship






Figure 2 Deconvolution of in vivo release profile (points and line correspond to raw and fitted data, respectively) and modelling of dissolution data (yellow points and line correspond to raw and fitted data, respectively)

Validation

- Only the non-linear approach met FDA validation criteria
- The % PE's of the non-linear approach are summarised in Table 1

Formulation Obs C_{max} Pred C_{max} C_{max} Obs AUC Pred AUC AUC %

Figure 3 Comparison of linear (top) and non-linear (bottom) IVIVC approaches. From left to right: fitted Levy plots, mean scaled in vitro profile (yellow points), predicted PK from established IVIVC

Conclusions

In this work, we have demonstrated a step-by-step approach for non-linear IVIVC using higher order polynomials.

The results showed that in this instance, when dissolution was much faster than absorption and the complexity of the release profile was high, a linear IVIVC was invalid and said to be inconclusive whereas a non-linear approach led to a valid IVIVC.

	(ng/mL)	(ng/mL)	PE	(ng*h/mL)	(ng*h/mL)	PE
1	59.7	56.8	4.73	26500	26400	0.200
2	54.2	52.0	4.08	24100	24200	0.400
3	64.5	59.5	7.75	29200	27700	5.24
MAPPE			5.52			1.95

Table 1 Observed (Obs) and predicted (Pred) Cmax and AUC following non-linear IVIVC and the associated percent error for each formulation.

References \

1. FDA Guidance for Industry Extended Release Oral Dosage Forms: Development, Evaluation and Application of In Vitro/In Vivo Correlations (1997)

Mauger DT & Chinchilli VM (1997) J Biopharm Stat. 1997;7:565–78. https://doi.org/10.1080/10543409708835207
Bhalan G, Timmins P, Greene DS, Marathe PH J Pharm Sci. 2001; 90:1176–85. https://doi.org/10.1002/jps.1071
Mohamed MEF., Trueman S, Othman AA et al (2019) AAPS J 21: 108. https://doi.org/10.1208/s12248-019-0378-y
Shen J, Burgess DJ (2015) J. Control. Release 219 644–651. https://doi.org/10.1016/j.jconrel.2015.09.052

6. Dickinson J, Patterson CM, Gajjar P, Dickinson H and Dickinson PA Development of a methodology to enable non-linear in vitro-in vivo correlation for complex long-acting injections *In preparation*



Pharmaceutical Development Services

E: enquiries@sedapds.com W: www.sedapds.com Seda Pharmaceutical Development Services® is the business name and registered trademark of SEDA Pharma Development Services Ltd, a company incorporated in England and Wales with registered number: 9442533 and registered office: 3 Castlebrook Close, Unsworth, Bury, Lancashire, UK, BL9 8JE. ©Copyright 2015-2017